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


Nuclear receptor

Binding assay
Cofactor recruitment assay
Reporter gene assay
Nuclear translocation assay


Ion channel / Transporter

Ion influx assay
Membrane potential
Electrophysiology
Substrate uptake
Binding

PPI (protein-protein interaction)

TR-FRET/Alpha screen assay
ELISA
NanoBit/BRET
Two-hybrid assay
Biophysical assay


J Med Chem **56** 9635–45 (2013)

Enzyme

Luminescence, Absorbance, Fluorescence, TR-FRET
Alphascreen, ELISA
Radiometric assay
Label-free assay (e.g. Rapidfire-MS)
Coupling assay
Global kinase panel

Nucleic acid

Biophysical assay (e.g. ASMS)
Fluorescence probe binding
FRET
Cell-based assay (Reporter gene, RT-qPCR)

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Enzyme

- Luminescence, Absorbance, Coupling, Fluorescence, TR-FRET, AlphaScreen
- Label-free assay (Rapidfire-MS)
- ELISA
- Radiometric assay

GPCR

- cAMP assay
- Ca²⁺ flux assay (FLIPR, FDSS)
- Reporter gene assay
- Arrestin/Internalization assays
- Binding assay

Phenotype

- High content imaging assay
- Reporter gene assay, qRT-PCR
- Cell growth etc.

Protein/Protein Interaction

- TR-FRET, AlphaScreen
- ELISA
- NanoBit, NanoBRET
- Two-hybrid assay
- Biophysical assay

Ion channel / Transporter

- Ion influx assay
- Membrane potential
- Electrophysiology
- Substrate uptake
- Binding

Biophysics

- SPR
- ASMS
- ITC
- Thermal Shift
- Crystal structure



Nucleic Acids

- Biophysical assay (e.g. ASMS)
- Fluorescence probe binding
- FRET
- Cell-based assay (Reporter gene, RT-qPCR)

Nuclear Receptor

- Binding assay
- Cofactor recruitment assay
- Reporter gene assay
- Nuclear translocation assay

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(Agilent)

Syncropatch 384 (Nanon)

QuantiStudio 12K Flex Realtime PCR
(ThermoFisher)

IN Cell Analyzer 6000
(Molecular Devices)

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


Compound No	Cluster ID	Primary assay IC50 (M)	Counter assay IC50 (M)	Primary assay (Graph)	Counter assay (Graph)	Chemical properties				Purity (%)	Notes
						MW	HBD	HBA	AlogP		
AXL1	1	8.0E-08	>1.0E-05			...	...	...	...	95.5	kinase inhibitor
AXL2	1	2.9E-07	>1.0E-05			...	...	...	...	90.5	kinase inhibitor
AXL3	1	4.9E-06	>1.0E-05			...	...	...	...	95.3	kinase inhibitor
AXL4	2	3.2E-08	>1.0E-05			...	...	...	...	85.4	promiscuous
AXL5	3	1.2E-06	>1.0E-05			...	...	...	...	97.2	
AXL6	3	2.5E-06	>1.0E-05			...	...	...	...	98.1	
AXL7	3	4.2E-06	>1.0E-05			...	...	...	...	95.4	

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

GPR39 positive allosteric

Primary screening


Library: >600,000 cpds at 3 μ M
hGPR39 cAMP assay (Gs) with EC₂₀ Zn²⁺

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

HTS cascade

Primary screening (ca. 370,000 compounds)

- Diversity pooled library 10 μ M, N=1
- Enzyme assay with Rapidfire-MS


Deconvolution assay (ca. 800 cpds)

- Positive compounds from primary screening 30 μ M, N=1


Dose response test

- 4-5 dose, N=2
- Selectivity test (LDHB/LDHA)
- Clustering
- Purity check
- Evaluation of related compounds

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HT-ASMSを用いたRNA binder探索

High-throughput Affinity Selection Mass Spectrometry for FMN riboswitch RNA binders

Diversity Library: 140,633 cpds
RNA focused Library: 6,400 cpds

FMN riboswitch RNA 5 μ M
Compound conc.: each 0.5 μ M
Assay format: 400 cpds/pool




↓
Reproducibility test

↓
Counter Screen using scramble RNA

↓
Confirmation assay

Diversity Library : **72** hits
RNA focused Library : **7** hits

↓
Kd determination

In house data

A β uptake assay in iCell[®] Microglia AD TREM2

HTS cascade

Primary screening (c.a. 4000 compounds)

- Biologically annotated compounds, > 3,000 cpds
- 3 μ M, N=1
- Phagocytosis assay and cytotoxicity (CellTiter-Glo[®] Luminescent Cell Viability Assay)

Reproducibility test (350 compounds)

- Positive compounds from primary screening
- 3 μ M, N=1
- Phagocytosis assay and cytotoxicity (CellTiter-Glo[®] Luminescent Cell Viability Assay)




Dose response test (24 compounds)

- Selected compounds from reproducibility test
- 6 dose, N=2
- TREM2 mutant and WT
- Phagocytosis assay and cytotoxicity (CellTiter-Glo[®] Luminescent Cell Viability Assay)

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

<Strategy>

<New finding>

Feedback activation of autophagy is a key resistance mechanism against SCD1 inhibitor-induced cell growth inhibition.

SCD1 inhibitor and autophagy inhibitor co-treatment

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